

MASSTEC MEDICAL CO.

GYNECOLOGICAL CATHETERS

- Intrauterine Insemination Catheter
- Ovum Aspiration Needle
- Embryo Transfer Catheter Set
- Endometrial Suction Curette
- Hysterometer
- Intrauterine Device Hook



Background

Masstec Medical Co. is a medical device manufacturing company, founded in 2013 as a subsidiary of Poursina Hakim Holding Co. Taking advantage of prominent physicians' guidance and deployment of expert engineers has led to the company gaining a deep understanding of **physicians' in-action challenges** and the **pain patients go through**.

1 🗌

Achievements:

• Achieving the technical knowledge and skills needed to design and manufacture ovum aspiration needles and embryo transfer catheters.

• Introducing the disposable models of some devices for the first time in Iran (products such as the company's disposable endometrial suction curette, hysterometer, and IUD hook).

• The company is among the Iranian market's leading manufacturers of disposable gynecology catheters.

Massiah's life-giving breath:

"The delight of healing" is Masstec's motto which concisely explains our organizational mission.



masstecmedical.com



Certificates and Permits









masstecmedical.com

4



— IUI Catheter

Intrauterine Insemination is considered one of the oldest and most practical methods for treating infertility. In this method, concentrated sperm is introduced into the uterus. For this purpose, a special catheter, called an IUI catheter, is used.

Models:


• Simple IUI Catheter and IUI Catheter with Mandrel

Specifications: Luer-lock hub **Stopper** • Dimensions: 17 cm x 8 Fr Ensure the syringe Measure the • EO sterilized connection insertion depth In compliance with 10993 SO ISO 11737 BS EN 1618 S S **Close-end** masstecmedical.com Atraumatic design



– IVF Needle – Ovum Aspiration

Oocyte (egg) retrieval or ovarian puncture is a procedure carried out in the IVF process during which the ovary is stimulated using specific drugs and then, using a special needle called the Ovum Aspiration Needle, and under the guidance of ultrasound imaging, oocytes are extracted from the ovary and stored in laboratory conditions.





IVF Catheter – Embryo Transfer

During the IVF procedure, the embryo is placed into the uterus using a special device called an Embryo Transfer Catheter. This device comprises two catheters, namely the guiding catheter and the transfer catheter.

Specifications:


- Guiding catheter: 17 cm × 6.6 Fr
- Transfer catheter: 24 cm \times 2.7 Fr





— Endometrial Curette – Simple Suction

Endometrial and uterine biopsy are generally performed on those suspected to have cancer in order to accurately diagnose the condition. Using an endometrial suction curette enables physicians to extract the required sample tissue for pathology tests by applying suction and gently scraping the endometrium, while causing minimal trauma and discomfort.





— Endometrial Curette – Syringe Suction

Specifications:

- Dimensions: 22 cm \times 9 Fr
- EO sterilized





— Hysterometer

A hysterometer is a device used for this purpose. Although previously reusable hysterometers were more widely used, nowadays, due to higher hygienic standards and to prevent the transmission of contagious diseases, disposable hysterometers are more in demand.

Specifications:

- Dimensions: 17 cm \times 12 Fr
- EO sterilized





d IUD Hook

This disposable device, which has two notches at its distal end, is used for extracting Intrauterine Devices (IUD) from the patient's uterus.





www.masstecmedical.com
info@masstecmedical.com
info@masstecmedical
+98 912 057 14 62